

EMERGENCY CARE · GP PRACTICES & CARE HOMES

Defibrillator (AED) weekly check log

A simple weekly record that your AED is ready to use: status light, pads, battery and the kit around it.

AED make and model _____ Serial number _____

Location _____ Pads expiry _____

Battery expiry _____ Responsible person _____

Week commencing	Status OK?	Pads sealed and in date?	Battery OK?	Kit complete?	Initials	Notes

— What "kit complete" means

- ☐ Adult pads (and child pads or a child key if you might need them)
- ☐ A spare set of pads
- ☐ Gloves, scissors, razor and a pocket mask
- ☐ A clean, dry towel

Top tip

After any use, replace the pads straight away and follow the manufacturer's instructions for downloading data and resetting the device.

Stock up in one order

Order by 4pm for next-day delivery. One account, one invoice.

[Which defib pads? →](#)

[Zoll defibrillators →](#)

This resource is general guidance to help you organise your practice. It does not replace national guidance, manufacturer instructions or your own policies, which you should always follow. Check that it reflects the current requirements that apply to your organisation. Produced by Medical Supermarket, a trading name of ProCure Health Ltd.